

## **Cryopreservation (droplet vitrification) of *Malus* using apical shoot tips derived from plants grown *in vitro***

**USDA-ARS National Laboratory for Genetic Resources Preservation, Fort Collins, CO**

### **Plant Material Culture Conditions**

*Malus in vitro* stock cultures are maintained in an actively growing state on multiplication medium in growth room at 25°C with a 16 h photoperiod (40  $\mu\text{M m}^{-2}\cdot\text{s}^{-1}$ ) and subcultured once every 4 weeks .

### **Shoot Tip Excision and Preculture**

Harvest apical shoot tips (1.5 mm) from 4-week-old *in vitro* stock cultures and plate on basal medium overnight at 25°C in darkness. Remove shoot tips from basal medium and place on preculture medium; preculture for 1 days at 25°C in darkness.

### **Cryopreservation Method**

Remove shoot tips from preculture medium and cryoprotect with full-strength PVS2 (Sakai et al. 1990) at room temperature for 20-40 minutes (40 minutes has been the most effective time for most species/cultivars). Place shoot tips into very thin layer of PVS2 on foil strips and plunge into liquid nitrogen.

### **Warming**

Warm foil strips in room temperature unloading solution and hold for 20 minutes. Plate shoot tips onto basal medium overnight at 25°C in darkness. Transfer shoot tips to fresh basal medium and culture for 7 days at 25°C in darkness and then expose to light at 25°C. Evaluate for survival/recovery after 8 weeks, transferring to fresh media if necessary.

### **References**

Li B, Feng C, Wang M, Hu L, Volk G, Wang Q. 2015. Recovery patterns, histological observations and genetic integrity in *Malus* shoot tips cryopreserved using droplet-vitrification and encapsulation-dehydration procedures. *Journal of Biotechnology* 214:182-191.

Sakai A, Kobayashi S, Oiyama I. 1990. Cryopreservation of nucellar cells of navel orange (*Citrus sinensis* Osb. Var. *Brasiliensis* Tanaka) by vitrification. *Plant Cell Rep.* 9(1):30-33. doi: 10.1007/BF00232130

Wang M, Chen L, Teixeira da Silva JA, Volk GM, Wang Q. 2018. Cryobiotechnology of apple (*Malus* spp.): development, progress and future prospects. *Plant Cell Reports* 37:689-709.

## Media and Solutions

- **Multiplication medium:** MS + 30 g L<sup>-1</sup> sucrose + 0.1 g L<sup>-1</sup> myo-inositol, 1 mg L<sup>-1</sup> 6-benzyl aminopurine (BAP), 0.3 mg L<sup>-1</sup> indole-3-butyric acid (IBA), 0.2 mg L<sup>-1</sup> gibberellic acid and 7.5 agar at pH 5.5 (pH 5.8 prior to autoclaving)
- **Basal medium:** MS + 30 g L<sup>-1</sup> sucrose + 0.25 mg L<sup>-1</sup> BAP + .01 mg L<sup>-1</sup> IBA + 2.6 g L<sup>-1</sup> gellan gum at pH 5.7 (pH 6.0 prior to autoclaving)
- **Preculture medium:** MS + 2M glycerol + 0.8M sucrose at pH 5.8 (pH 6.4 prior to autoclaving)
- **PVS2 solution:** 1/2 MS + 30% (w/v) glycerol + 15% (w/v) ethylene glycol + 15% (w/v) dimethyl sulfoxide + 0.4 M sucrose at pH 5.8 (filter sterilized solution)
- **Unloading solution:** MS + 1.2M sucrose at pH 5.7 (pH 6.6 prior to autoclaving)
- **Other parameters:** wrap all cultures with two layers of plastic wrap (not Parafilm); after pouring plates allow to cool for thirty minutes before covering/packaging